

# YOUR MANAGED EMERGENCY SOLUTION SERVICE

**Assist me: Is activated in 2 ways. (a)** By pressing and holding and releasing the SOS button for 2 seconds.  
**(b)** Should the wearer slip/ fall the device will self activate.  
**Is** your connection to medical help when you need it most.  
**Works** anywhere on a national basis using cellular towers.

The emergency control center which is manned 24 hours a day, 7 days a week, 365 days a year will receive your GSM GPS location, you will then receive a call directly to your Assist me panic button.

2 SMS's will be sent to 2 watchers of the unit. Sending them a notification that the panic button has been pressed, it will also send them the device location, unit contact number and access to google maps to navigate the watcher to the device location.

The medical dispatcher will enquire the nature of the emergency and then dispatch the relevant emergency service to the clients location

<b>Assist Me Panic Button.</b>	
Once off fee: R 1800.00 (Excluding VAT)	
Monthly recurring fee: R 149.00 (Excluding VAT)	
Access to our 24hour medical call center*	
Single button or fall activation	
Two way voice communication with the unit itself	
Emergency SMS notification sent to 2 watchers	
Location of the pressed button sent to watchers	
Battery life of 5-6 days	
Low battery SMS notification	
Splash proof (IP65)	
Charging cradle	
7 day non- satisfactory money back guarantee	
12 month guarantee	

\*Ts & Cs Apply